

CartiBlaster™ Rasps

Powered Joint Prep Instruments

Prep joints in a
fraction of the time.
And have a blast doing it!



Experience the speed and control of CartiBlaster™ powered joint prep
and never go back to manual prep again.

TREACE®

CartiBlaster™ Rasps

Powered Joint Prep Instruments

CartiBlaster Fan™

Fast and efficient

with double-sided design for dual-joint prep

Open window

for efficient tissue evacuation to prevent clogging

Anatomic curvature

designed to maintain the natural contours of the mid and hindfoot joints

Applications: Subtalar, TN, CC, NC, TMT



CartiBlaster Edge™

Precision cutting edge

for controlled joint debridement and bone contouring

Curved shaft

allows for access to challenging curved joint surfaces

Highly versatile design

for application across multiple joints and bony surfaces

Applications: Lateral side of talar head, concavity of navicular, medial aspect of subtalar, MTP

Sterile-packaged and sharp for every case.

Ordering Information

SN57	CartiBlaster Fan™ - SM
SN58	CartiBlaster Fan™ - ZMS
SN59	CartiBlaster Edge™ - SM
SN60	CartiBlaster Edge™ - ZMS

Before use of the instruments, the surgeon should refer to the appropriate instructions for use for complete warnings, precautions, indications, contraindications, and adverse events. Risks include, but are not limited to: infection, pain, discomfort, nerve or soft tissue damage, and necrosis of tissue or inadequate healing. If any of these occur, additional treatments may be needed. Additional information about risks, warnings, and instructions is available at [Treace.com/labeling](https://www.treace.com/labeling).

**Order CartiBlaster™
for your next fusion.**



TREACE
A Step Ahead.™