TreaceForward™

BarrelReamer™

SpeedMTP® Precision Joint Prep Instruments

Streamline Your MTP Joint Prep.



Designed to efficiently remove all metatarsal exostosis, in a single pass.

TREACE®

BarrelReamer™

SpeedMTP® Precision Joint Prep Instruments


Expansive Viewing Window

Offers clarity for accurate contouring and alignment



Inward-Curving Blades

Minimizes tissue interference and enhances precision




Consistently Sharp


Sterile-packed and compatibile with Treace® 1st MTP Reamers.

Deep barrel depth allows for thorough and efficient removal of all osteophytes down to the 1st metatarsal shaft.

BEFORE



AFTER



Ordering Information:

Part#	Size
SN50	18mm BarrelReamer™ Instrument
SN51	21mm BarrelReamer™ Instrument
SN52	23mm BarrelReamer™ Instrument

For more information, go to **Treace.com**

TREACE® Medical Concepts, Inc.

Before use of the instruments, the surgeon should refer to the appropriate instructions for use for complete warnings, precautions, indications, contraindications, and adverse events. Risks include, but are not limited to: infection, pain, discomfort, nerve or soft tissue damage, and necrosis of tissue or inadequate healing. If any of these occur, additional treatments may be needed. Additional information about risks, warnings, and instructions is available at Lapiplasty.com/surgeons/labeling.

Pat. treace.com/patents. | ©2026 Treace Medical Concepts, Inc. All rights reserved. M4783A